

Aspiration Needles

for Bronchoscopy & Gastroscopy - Single Use

Aspiration Needles

- Coaxial Luer-Lock for Vacuum Syringe
Easy Aspiration of Cells
- Safety Clip and Lock Mechanism
Secure Needle Advancement and Projection
- 15° Beveled Needle with Side Hole
Atraumatic Piercing Capability
Maximum Yield of Aspirates



Aspiration Needles

for Bronchoscopy & Gastroscopy - Single Use

Bronchoscopy

Aspiration Needle – with Metal Spiral and including 10 ml Vacuum Syringe with Stopcock

Catalog No.	Needle Diameter	Needle Length	Metal Spiral Outer Diameter	Working Length	Minimum Accessory Channel
GAN-01-19-100	21 Gauge / 0.8 mm	12 mm	1.9 mm	100 cm	2.1 mm

Gastroscopy

Aspiration Needle – with Plastic Tube and including 10 ml Vacuum Syringe with Stopcock

Catalog No.	Needle Diameter	Needle Length	Sheath Outer Diameter	Working Length	Minimum Accessory Channel
GAN-02-19-180	21 Gauge / 0.8 mm	12 mm	1.9 mm	180 cm	2.1 mm

Distributed by:

Manufactured by: