

Retrieval Devices


Single Use

Polyp & Foreign Body Retrievers


- Strong & Wide Opening Graspers
Large variety of applications
- Ergonomically Shaped Handle
User-friendly and safe

Foreign Body Forceps

- Coated Metal Spiral with Positioning Bands
Friction minimized passage through working channel
- Convenient evaluation of forceps position for a controlled & safe catheter advancement and removal*



Polyp & Foreign Body Retrievers



Foreign Body Retriever – Grasper with Hooks

Catalog No.	Grasper shape	Catheter outer diameter	Working length	Minimum working channel
GFB-01-25-180	4-Prong Grasper	2.5 mm	180 cm	2.8 mm
GFB-01-25-230	4-Prong Grasper	2.5 mm	230 cm	2.8 mm

Polyp Retriever – Grasper with Eyelets

Catalog No.	Grasper shape	Catheter outer diameter	Working length	Minimum working channel
GFB-06-25-180	3-Prong Grasper	2.5 mm	180 cm	2.8 mm
GFB-06-25-230	3-Prong Grasper	2.5 mm	230 cm	2.8 mm
GFB-02-25-180	4-Prong Grasper	2.5 mm	180 cm	2.8 mm
GFB-02-25-230	4-Prong Grasper	2.5 mm	230 cm	2.8 mm

Foreign Body Retriever – Grasper with Forked Jaws

Catalog No.	Grasper shape	Catheter outer diameter	Working length	Minimum working channel
GFB-03-25-180	Forked Jaws 2:2	2.5 mm	180 cm	2.8 mm
GFB-03-25-230	Forked Jaws 2:2	2.5 mm	230 cm	2.8 mm

Foreign Body Basket

Catalog No.	Grasper shape	Basket diameter	Catheter outer diameter	Working length	Minimum working channel
GFB-07-25-425	4-Wire Basket	25 mm	2.5 mm	230 cm	2.8 mm
GFB-07-25-625	6-Wire Basket	25 mm	2.5 mm	230 cm	2.8 mm

Foreign Body Forceps

Foreign Body Forceps



Rat Tooth



Alligator

Foreign Body Forceps – with Coated Metal Spiral


Catalog No.	Cup shape	Catheter outer diameter	Working length	Minimum working channel
GFB-09-23-120	Rat Tooth	2.3 mm	120 cm	2.5 mm
GFB-09-23-180	Rat Tooth	2.3 mm	180 cm	2.5 mm
GFB-09-23-230	Rat Tooth	2.3 mm	230 cm	2.5 mm

Foreign Body Forceps – with Coated Metal Spiral

Catalog No.	Cup shape	Catheter outer diameter	Working length	Minimum working channel
GFB-10-23-180	Alligator	2.3 mm	180 cm	2.5 mm
GFB-10-23-230	Alligator	2.3 mm	230 cm	2.5 mm

Retrieval Devices

Single Use



Distributed by:

L-GFB-SU-00-ENG Ver. 5.0 (2016-02)